



1. **Instep Strap** assists hindfoot control. Maintains brace's positive closure, and should be worn snugly. For least bulk, order velcro only. For greater durability, use webbing-backed velcro.
2. **Forefoot Strap** helps secure the foot and insures good closure of the brace. It is optional, but recommended, especially for larger braces.
3. **Toe Abduction Strap** exerts tremendous influence on the toes and forefoot. Especially useful in severe pronation. Does not allow the great toe to abduct into eversion of the toes and forefoot.
4. **Anterior Stop (or Pre-Tibial) Strap** helps proprioception. Use with positioning braces or braces used in standing frame and prone standers. PLEASE NOTE: This strap is not strong enough to block the tibia from flexing forward. DO NOT USE TO RESIST A CROUCHED GAIT (Use Rear Entry Floor Reaction DAFO #1).
5. **Mild Plantarflexion Stop Strap** is useful when working with a child who requires augmented plantarflexion control during certain activities, but is beginning to develop some graded control. Therapist can engage or disengage strat during therapy quickly, without removing the splint.


TIP: Toe abduction strap may cause some patients discomfort by pulling sock down on ends of toes and into the space between the 1st and 2nd toes. To alleviate this problem, sew a sharp "V" in the sock between the toes. Then split the "V" with scissors to create a foot "mitten" with space for the great toe.



FREE PLANTARFLEXION
DAFO #4



RESIST PLANTARFLEXION
DAFO #3 ½



STOP PLANTARFLEXION
DAFO #3