

## Cascade Dafo, Inc.

1360 Sunset Ave, Ferndale, WA 98248 ph 800.848.7332 | intl +13605439306 fax 855.543.0092 | www.cascadedafo.com



PF resist, DF free, Softy liner

_	Last name:		Construction   Features   Options
ien	First name:		
Patient	Birth date: Bilateral Left Right		MEDIAL (Left) LATERAL (Left)
	Date cast:		
er	Last Name:	\ <del>\ \ \ \</del>	Proximal Strap
Practitioner	First Name:	$\Box$	Liner
i i	Title/Credentials:		Instep & Forefoot Straps
Pra	Email: Phone:		
50	Customer/Business name:		
	Street address:		
Billing		NOTE:	If you don't choose an option, you will receive the <b>Standard</b> .
m	City: State: Zip:	Liner:	Softy foam (white only) Only option available
	PO# / UCAN#:		Add extra navicular padding
	☐ Shipping info is the same as Billing info. —OR—	Ctrono	(boney pronators)
<b>.</b>	Shipping contact name:	Straps:	Standard (see drawing) Add toe abduction strap
Shipping	Facility name:	Strap Color:	White (Standard) Other:
hip	Street address:		NOTE: The posterior strap is available in white dacron only.
S		Transfer Pattern:	No Transfer (Standard)
	City: State: Zip:	rattern.	
	Finished Brace Angles		Pattern:  NOTE: Outer frame only; additional cost per brace.
ANKLE ALIGNMENT (Dorsiflexion-Plantarflexion)  Toe Rise and			
Correct to 3-4° DF Correct to OD not correct (cast alignment OK) Do not correct (cast alignment OK) Color: White (Standard) Other:			
HINDFOOT ALIGNMENT			
Correct to vertical (if misaligned)  Do not correct  Toe Shelf			
FOREFOOT ALIGNMENT NOTE: Drawings show finished orthosis Inner Liner: Medial			
	Choose forefoot alignment. Write posting height if needed- in. or mm.		no containment containment Standard
RIGHT	RICHT ROOM		AND / OR  Lateral containment
	Alana Mantal Mantal Mana		
_	/algus Varus Neutral Neutral Varus Valgus		
ᆜ			
	Do not correct - keep as cast.   Do not correct - keep as cast.		
Bottom Stabilization		Special Instructions	
	None- Standard		
	Heel -OR- Midfoot -OR- Both -OR-		
	Entire bottom stabilized with foam sole -OR-		
	Entire bottom with non-skid cover <b>-OR-</b>		
	Entire bottom stabilized with both foam sole <b>and</b> non-skid cover		
N	OTE: Varus or valgus forefoot alignments will receive stabilization on bottom		
NI/	of brace to support posted (raised) region		
M	OTE: Neutral forefoot alignments will not see foam on toe shelf		
		Rush ord	er (adds \$20)