





WHAT is it?

The JumpStart Softback is a prefabricated AFO providing comfortable control of mild to moderate toe walking and knee hyperextension. The Softback's unique posterior strap and pad combination work together to strongly resist plantarflexion while encouraging increased dorsiflexion.




Soft yet durable posterior padding sewn to flexible Dacron strap strongly resists excess plantarflexion and encourages free dorsiflexion.

Dacron strap attaches to outer frame for durability.



Full-wrap inner liner provides alignment control of the heel, midfoot, and forefoot.




Inner liner tucks behind the padding for added comfort.


Solid outer foundation provides additional support and heel stabilization.

WHO is it for?

The JumpStart Softback is designed to correct mild to moderate excess plantarflexion. It's an ideal solution for early intervention patients with fully-correctable foot position and some voluntary control.



Toe walking correction



Hyperextension correction

HOW does it help?

The Softback's comfortable pad and strap combination work together to strongly resist plantarflexion beyond 2 degrees, while encouraging increased dorsiflexion. This forgiving resistance allows easier access to crawling and floor-to-standing activities.

Sizing ?

This prefabricated brace can be sized by measurement using the Fast Fit sizing jig.

Widths: Wide and Narrow

Sizes: 4.00–9.00 (0.25-inch increments)

To order: contact Customer Service at 800.848.7332 or customerservice@dafo.com.