

The DAFO 8 Softy uses a soft foam liner to provide comfort and good biomechanical, long-term positioning for non-ambulatory patients. The orthosis provides full wraparound foot support for excellent hindfoot and forefoot control, and its proximal upright blocks both plantarflexion and dorsiflexion. The DAFO 8 Softy is suitable for supported standing transfers, yet is designed to provide superior comfort in non-ambulatory situations. The orthosis is ideal for patients who require a strong amount of control to maintain a proper foot and ankle position.



Key Features

Custom-fabricated from a cast or scan

Precise positioning

for non-ambulatory patients; suitable for supported standing and transfers

Defined plantar surface contours

Softy Liner

made of foam provides more comfortable control of feet that are difficult to align or have sensitive bony prominences

Designed for patients who:

- Are non-ambulatory.
- Have strong fixing patterns of the foot and ankle.
- Need wraparound foot control in order to maintain proper positioning
- Need very precise but compliant control.
- Require extra padding for feet with sensitive bony prominences.