


# Fast Fit® Chipmunk®



CASCADE  
daFO®

*Prefabricated maximum support shoe insert*

Designed to control moderate low tone pronation with an optimal amount of support and comfort, the Chipmunk has an integrated two-part system combining internal support from molded foam contours with the stability of a flat, external base and deep heel cup for maximum correction of foot position.



## Key Features

### High medial trimline

covers the navicular to protect from rubbing.

### Bottom plastic base

with diagonal forefoot trimline to facilitate natural push-off. Full heel cup for stabilization.

### Built-in arch padding

gives additional internal support.

### Plantar surface

supports good foot alignment with soft, resilient contours.

### Foam liner

comfortable and very flexible.

### Flattened toe shelf


leaves room for toes within the shoes, and can be trimmed to fine-tune the fit. Allows for addition of toe rise pad when desired.

### Prefabricated

sized by measurement using the Fast Fit sizing jig.

**Widths:** Wide and Narrow

**Sizes:** 4.00 – 12.25 inches (0.25 inch increments)



*Moderate pronation, unstable foot position.*

*Improved foot alignment and steadier gait.*

## Designed for patients who:

- Exhibit moderate pronation and have a fully correctable foot position. The patient's medial arch is reduced, forefoot is moderately abducted, and heel is everted.
- May benefit from the sensory feedback that a contoured shoe insert can provide.



**No casting needed!**

Chipmunk product sheet rev.07 (Sep 2024)

Helping Kids lead healthier, happier lives®

Cascade Dafo, Inc.

1360 Sunset Avenue  
Ferndale, WA 98248

ph: 800.848.7332  
fax: 855.543.0092

intl: +1 360 543 9306  
[www.cascadedafco.com](http://www.cascadedafco.com)